DESIGNED  
**USA**  
PRODUCED

# Simplisse™

Manual  
Breastfeeding  
Companion™

INSTRUCTION BOOKLET

Manual Breastfeeding  
Companion  
*Instructions For Mom*

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
## THE SIMPLISSE STORY

Now, giving your newborn the benefits of nutritional breastmilk doesn't have to be a painful experience. After years of breastfeeding research and observation, a team of leading lactation consultants partnered with moms to create a milk collection system that simulates the comfort and benefits of baby at the breast. Simplisse Breastfeeding Companions feature a unique BabyFace™ Flexible Breastcup that authentically emulates baby's mouth movements. It conforms to mom's breast just like baby, and it alleviates the discomfort and friction often associated with traditional suck-and-release breast pumps. In addition, Continuous Comfort™ Valves mimic baby's feeding pressure without painful tugging or pinching of mom's nipple. And when mom is comfortable and happy, breastmilk quality is at its highest. Ultimate comfort for mom and optimal milk flow for baby—these are the breakthrough benefits of the Simplisse Breastfeeding Companion.

**Simplisse: Naturally Better. By Design.™**

### TABLE OF CONTENTS:

Parts Diagram	3
Disassembly	4-5
Cleaning	5
Assembly	6-8
Expressing Breastmilk	8-9
Storing Breastmilk	10
Preparing Breastmilk	10
Tips	11
Troubleshooting	12

# PARTS DIAGRAM

- A. Stand
- B. Collection Bottle
- C. Valve
- D. Connector
- E. Breastcup
- F. Stem
- G. Large Insert
- H. Small Insert
- I. Cap
- J. Handle

## DISASSEMBLY

Follow these steps to disassemble your pump for cleaning:

1

First loosen and remove the Collection Bottle from the Connector.

2

Then remove the white Valve from the underside of the Connector.

3

To remove the Breastcup, gently pull on the perimeter tabs to separate it from the Connector.

4

To remove the Handle, first push the Stem loose from the slot on the Handle with your finger.

5

Then gently twist the Handle to separate from the post on the Connector.

6

To remove the Cap, simply pull up on the Cap tab at the back of the pump. This will separate from the Connector.

7

With the Cap in hand, pull the Stem free by pulling downward.

8

Pull the Large Insert from the Cap.

9

Finally, pull the Small Insert from the Cap.

**Note:** This pump has two small valves in the Connector. Do not remove these valves. These are the Continuous Comfort™ Valves which regulate the pump's pressure so there is no high suction. These are permanent parts.

## CLEANING

Before your first use and after all other uses, it is important to wash or sterilize your pump:

Make sure all pump parts have been disassembled for cleaning or sterilization.

All pump parts can be hand-washed, steam sterilized, boil sterilized, or washed in a dishwasher.

**Note:** The white Valve and Breastcup are delicate. Please take extra care in cleaning them. It is recommended that you hand wash these parts. Do not use any rough brushes to clean these parts.

**TO HAND WASH:** Wash all parts in warm, soapy water, rinse well and air dry.

**TO STEAM STERILIZE:** Use a microwave steam sterilizer or steam sterilizer bags. Be sure to follow the instructions that come with your sterilizer or steam bags.

**TO BOIL STERILIZE:** Boil all parts for 5 minutes, carefully remove and air dry.

**IN DISHWASHER:** All pump parts are top-rack dishwasher safe. Completely air dry parts and let them cool.

## ASSEMBLY

After cleaning, follow these steps to assemble your pump. Make sure all parts are completely dry from cleaning before any assembly.

1

Begin by placing the white Valve into the bottom of the Connector with the wedge side pointing down.

2

Then secure the Collection Bottle to the Connector. Tightening the Collection Bottle firmly will seal the Valve but do not over tighten.

3

Install the Breastcup by pushing at the edges to secure to the Connector. Make sure the Breastcup is sealed all around.

4

To prepare the Cap, first take the Small Insert and push it onto the small ring on the underside of the Cap.

5

Take the Stem and insert it through the Large Insert until it bottoms out.

6

The flat, round portion of the Stem should be flush with bottom of the Large Insert.

7

To complete the Cap, take the Stem and Large Insert assembly and push the head of the Stem through the Small Insert and through the top of the Cap.

8

Make sure the Large Insert is fixed to the large ring on the underside of the Cap. Aligning the Large Insert over the two small posts will orient it properly.

9

Take the Stem and pull up gently until the Small Insert snaps itself to the middle portion of the Stem.

10

Install the Cap assembly to the Connector assembly by pivoting the front latch to the Connector.

11

Then push down on the back portion to snap the tab on the Connector.

12

Install the Handle by catching the Stem in the slot of the Handle.

13

Then firmly push the Handle down to snap onto the post of the Connector.

14

**Note:** Make sure all parts are assembled correctly and snugly fit. The pump must be air tight to function properly.

## EXPRESSING BREASTMILK

Our pump is a very different design and does not work like traditional suck-and-release breast pumps.

### BabyFace™ Flexible Breastcup

We designed this pump with a unique breastcup that is very soft and flexible. It will conform to and move with your breast as you express. Our moms have told us that the flexible breastcup makes for a more comfortable pumping session. There is no friction or tugging, which can cause trauma to your breast, especially over time.

### Continuous Comfort™ Valves

Our pump was also designed with two small valves that constantly regulate the pressure of the pump so the suction levels are similar to that of baby at the breast. This means you don't have to worry about high suction and painful pumping.

Once your pump is cleaned and assembled, you are ready to begin pumping.

- 1 Make sure your hands and breasts are clean.

2

Position the Breastcup, nipple centered, to your breast just firm enough to ensure a complete seal with no air leaks.

3

Begin pumping by squeezing the Handle.

**Note:** You may hear an airy noise as you pump the Handle. This comes from the Continuous Comfort Valves and is completely normal.

Moms often achieve Let-Down with short, quick squeezes as these are similar to the fast, non-nutritive sucking motion of baby during the first phase of breastfeeding. This should initiate milk flow which is known as Let-Down.

### **Let-Down is the biological trigger to release breastmilk and it's the first step for a successful pumping session.**

Once you have Let-Down, continue to pump at a rate that is comfortable for you. Longer, slower squeezes often yield the most milk flow as these are similar to the actual nursing phase of breastfeeding.

Breastmilk should begin to flow and collect in the bottle. Even if you do not experience milk flow immediately, continue pumping for a few minutes. Your first few times expressing may not yield any milk. This is normal. Never exceed 5 minutes of pumping without milk flow.

Always clean or sterilize your pump after every pumping session.

**Warning:** Pumping should not cause pain or discomfort. If you experience discomfort during pumping, contact a Lactation Consultant or your healthcare provider.

4

5

Press side of breast to break suction.

## STORING & HANDLING BREASTMILK

After expressing your breastmilk, it is important to store and handle it properly.

- ⦿ Make sure your hands and storage containers are clean to avoid contamination of your expressed breastmilk.
- ⦿ Use Breastmilk Collection Bottles or Breastmilk Storage Bags to store your breastmilk. These can be used in the refrigerator or freezer.
- ⦿ Be sure your storage bottle or bag is sealed tightly to ensure freshness of the breastmilk.
- ⦿ Date code your bottles or bags if you do not intend to use your collected breastmilk immediately.
- ⦿ For extended storage, freeze your breastmilk.

**Warning:** Do not add fresh breastmilk to already frozen breastmilk in a storage container. Do not save breastmilk from a used bottle for another feeding.

## THAWING & PREPARING BREASTMILK

After storing your breastmilk, follow these steps to prepare it for baby:

**Warning:** Do not microwave breastmilk. Do not microwave Collection Bottles or Storage Bags containing breastmilk. Microwave ovens can create uneven heating which could scald baby and deteriorate the nutrients of breastmilk.

- ⦿ To thaw breastmilk, place the sealed Collection Bottle or Storage Bag in refrigerator or run under cool tap water.
- ⦿ To warm breastmilk, place the sealed Collection Bottle or Storage Bag in warm tap water.
- ⦿ Always check breastmilk temperature before feeding baby.
- ⦿ Use thawed breastmilk immediately or refrigerate up to 24 hours.

**Note:** Do not thaw breastmilk in very hot or boiling water. Once thawed, do not refreeze or refrigerate breastmilk again. Do not allow thawed breastmilk to sit at room temperature.

## Storage Duration of Fresh Breastmilk (For Use with Healthy Full Term Infants)

LOCATION	TEMPERATURE	DURATION
Countertop, table	Room temperature (up to 77°F or 25°C) <i>Containers should be covered and kept as cool as possible; covering the container with a cool towel may keep milk cooler.</i>	<b>6-8 hours</b>
Insulated cooler bag	5-39°F or -15-4°C <i>Keep ice packs in contact with milk containers at all times, limit opening cooler bag.</i>	<b>24 hours</b>
Refrigerator	39°F or 4°C <i>Store milk in the back of the main body of the refrigerator.</i>	<b>5 days</b>
Freezer compartment of a refrigerator	5°F or -15°C	<b>2 weeks</b>
Freezer compartment of refrigerator with separate doors	0°F or -18°C	<b>3-6 months</b>
Chest or upright deep freezer	-4°F or -20°C <i>Store milk toward the back of the freezer, where temperature is most constant. Milk stored for longer durations in the ranges listed is safe, but some of the lipids in the milk undergo degradation resulting in lower quality.</i>	<b>6-12 months</b>

Reference: Academy of Breastfeeding Medicine. (2004) Clinical Protocol Number #8: Human Milk Storage Information for Home Use for Healthy Full Term Infants. Princeton Junction, New Jersey: Academy of Breastfeeding Medicine.

## TIPS

Breastfeeding and pumping are unique experiences for nursing moms. Please review the tips below before pumping.

- ⊙ Moms are more likely to successfully pump if they are relaxed and thinking about baby. Avoid pumping if you are rushed.
- ⊙ Efficient pumping usually takes patience and practice for moms. Moms that pump frequently become accustomed to it and are more likely to have faster Let-Down and more milk flow.
- ⊙ All moms are different. Some experience Let-Down quickly and some moms do not. Some express on their first try and some moms do not.
- ⊙ Do not be discouraged if you are unable to express milk on your first try. Simply take a break and try again later.
- ⊙ Do not be discouraged if you express only a small amount of breastmilk. This is common especially on the first few attempts.
- ⊙ Pumping consistently, in addition to breastfeeding, can help sustain your milk supply over time.
- ⊙ A warm compress or a warm bath/shower can help encourage milk flow.
- ⊙ Most important of all—pumping should be a comfortable experience. If you feel pain or discomfort during pumping, contact a Lactation Consultant or your healthcare provider.

## TROUBLESHOOTING

- ⦿ For your pump to work as designed, all parts must be properly assembled. Make sure all parts are snugly fit so you have an air tight pump.
- ⦿ When applying to your breast, make sure the Breastcup is centered around your nipple. Hold the pump firm enough to your breast for a complete seal all around, but not too firm as this could inhibit milk flow.
- ⦿ When the Handle is squeezed, make sure that you feel a light suction on your breast.
- ⦿ If you hear an airy noise as you pump the Handle, this comes from the Continuous Comfort Valves™ and is completely normal.
- ⦿ If you require personal assistance about your pump, or have lost or broken parts, please call our Simplisse Customer Service Department at 877-962-2525.

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**Simplisse: Naturally Better. By Design.<sup>™</sup>**

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© 2009 Simplisse, Inc.