E

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GB

Instruction manual for BMI calculator seca model 491

The seca calculator model 491 will be used for the accurate determination of the Body Mass Index BMI. In addition to the BMI calculation the calculator offer the 4 basic arithmetics. The operation is equal to common calculators. Please care for sufficient illumination of the solar cell to avoid low contrast or a black display. For BMI calculation 8 different procedures can be used.

Activate the calculator **ON/C**

Procedure 1 : weight in kg / body height in m e.g. 80.5 kg / 1.82 m

Keyboard entry: **KG** **80.5** **M** **1.82** **BMI** = 24.3

Procedure 2 : weight in kg / body height in feet e.g. 80.5 kg / 5.97 ft

Keyboard entry: **KG** **80.5** **FT** **5.97** **BMI** = 24.3

Procedure 3 : weight in kg / body height in inches e.g. 80.5 kg / 71.65 in

Keyboard entry: **KG** **80.5** **IN** **71.65** **BMI** = 24.3

Procedure 4 : weight in kg / body height in feet / inches e.g. 80.5 kg / 5 ft 11.64 in

Keyboard entry: **KG** **80.5** **FT** **5** **IN** **11.64** **BMI** = 24.3

Procedure 5 : weight in lbs / body height in m e.g. 177.5 lbs / 1.82 m

Keyboard entry: **LBS** **177.5** **M** **1.82** **BMI** = 24.3

Procedure 6 : weight in lbs / body height in ft e.g. 177.5 lbs / 5.97 ft

Keyboard entry: **LBS** **177.5** **FT** **5.97** **BMI** = 24.3

Procedure 7 : weight in lbs / body height in inches e.g. 177.5 lbs / 71.65 in

Keyboard entry: **LBS** **177.5** **IN** **71.65** **BMI** = 24.3

Procedure 8 : weight in lbs / body height in ft and inches e.g. 177.5 lbs / 5 ft und 11.64 in

Keyboard entry: **LBS** **177.5** **FT** **5** **IN** **11.64** **BMI** = 24.3

***Note:**

-After BMI has been displayed input memory is cleared for a new input.

-If the unit can no longer be switched on, the unit's software has "crashed". If this happens use a sharp object such as a pencil to press the reset knob on the back of the unit.

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GB

Body Mass Index (BMI)
Body Mass Index puts height and weight in relationship to one another:

Body Mass Index (BMI) = weight in kg / (height in m)²

On the seca 491, even if you cannot enter metric measurements, the device will automatically convert to metric units.

Note: Body Mass Index applies only to adults: these values are not relevant to children and young people.

Evaluating Body Mass Index:

Compare the value determined with the categories quoted below which correspond to those used by WHO, 2000.

- BMI below 18.5:

Warning, this patient weighs too little. There could be a tendency to anorexia. An increase in weight is recommended to improve well-being and performance. If in doubt, consult a specialist.

- BMI between 18.5 and 24.9:

The patient is of a normal weight. He can stay as he is.

- BMI between 25 and 30 (pre-obese):

The patient is slightly to moderately overweight. He should cut his weight if he is already ill (e.g. diabetes, hypertension, gout, fat metabolism disorders).

- BMI over 30:

Weight reduction is essential. Metabolism, cardiovascular system and bones are all under strain. A consistent diet, plenty of movement and behavioural training are all recommended. If in doubt, consult a specialist.

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